



Adam's OIML weights are available in several classes.

E2 weights offer superb accuracy for calibrating analytical and precision balances.

Features

- Stainless steel weights for durability and accuracy
- Supplied in a plastic case for safe storage.

Individual stainless steel weights

Model	Class	Density	Susceptibility	Tolerance	Value
E2 50g Calibration Weight	OIML E2	7.96 g/cm ³	0.01 - 0.005	0.10 mg	50g
E2 100g Calibration Weight	OIML E2	7.96 g/cm ³	0.01 - 0.005	0.16 mg	100g
E2 200g Calibration Weight	OIML E2	7.96 g/cm ³	0.01 - 0.005	0.30 mg	200g

Accessories

UKAS calibration certificate for E2 weights

About ADAM:

For more than 45 years, professionals worldwide have relied on Adam Equipment for an extensive selection of dependable and affordable weighing equipment. Adam provides the right balance of speed, performance and value.

Supplied by