



Adam's OIML weights are available in several classes.

E1 weights are designed to be used as a reference for other calibration weights. They can also be used to calibrate semi-micro balances.

Features

- Stainless steel weights for durability and accuracy
- Supplied in a plastic case for safe storage.

Individual stainless steel weights

Model	Class	Density	Susceptibility	Tolerance	Value
E1 5g Calibration Weight	OIML E1	7.9 g/cm ³	< 0.5	0.016 mg	5g
E1 10g Calibration Weight	OIML E1	7.9 g/cm ³	< 0.5	0.020 mg	10g
E1 20g Calibration Weight	OIML E1	7.9 g/cm ³	< 0.5	0.025 mg	20g
E1 50g Calibration Weight	OIML E1	7.9 g/cm ³	< 0.5	0.030 mg	50g

About ADAM:

For more than 45 years, professionals worldwide have relied on Adam Equipment for an extensive selection of dependable and affordable weighing equipment. Adam provides the right balance of speed, performance and value.

Supplied by